

**European survey of C4d immunostaining protocols on frozen sections (total number of European centres: 14)**

		N. of European centres
<b>C4d antibodies</b>	Quidel (mc)	13
	Ventana 760-2687 (pc)	1
<b>Dilution</b>		
	Quidel 1/10 -1/500 (most used 1/100)	12
	Ventana 760-2687 : 1/25	1
<b>Incubation time</b>	1 hour	3
	5 – 45 minutes	7
	Overnight	1
	2 hours	2
<b>Incubation temperature</b>		
	Room temperature	8
	37°	2
	4°	2
<b>Fixation</b>	frozen	14